

Berikut contoh langkah-langkah menganalisis pola XRD dengan nama file “bajaST37.RD”. Ekstensi RD, merupakan file pola XRD dari mesin PANalytical. Tiap mesin XRD mempunyai tipe file tersendiri.

NOTE:

Untuk mengkonversi file pola XRD ke tipe yang lainnya dapat menggunakan aplikasi POWDLL.

Berikut langkah-langkah untuk menganalisis pola XRD menggunakan metode Rietveld dengan aplikasi X'PERT HIGHSCORE PLUS.

## 1. SEARCH-MATCH

Pada langkah ini, akan dilakukan pencocokan pola XRD dari sampel dengan pola XRD dari database. Database pola XRD ada beberapa, yaitu JCPDS/ICDD, COD dan lain-lain. Saat pencocokan pola XRD yang menjadi acuan adalah sudut/ d spacing puncak dengan intensitas tertinggi dan komposisi yang sebelumnya diperoleh dari XRF atau EDS (bila sampel batuan/ bahan yang tidak dikenal)

Dari database diperoleh data sistem kristal, space group, parameter kisi, posisi atom dan hkl dari puncak-puncak pola difraksi. Database ICDD PDF2 tidak menyediakan data posisi atom, sehingga diperlukan sumber lain, sedangkan COD sudah disertakan namun jumlah kartu pola difraksi lebih sedikit dibandingkan ICDD. Database ICDD yang menyediakan data posisi atom adalah PDF4, namun lisensinya hanya berlaku satu tahun saja, yang berbeda dengan PDF2.



Program yang digunakan adalah X'PERT HIGHSCORE PLUS berikut langkah-langkahnya.


- a. Menjalankan program XHP




Klik dua kali pada icon tersebut pada desktop untuk menjalankan XHP



- b. Membuka file pola XRD






c. Menentukan puncak-puncak pola XRD sampel.







Structure Plot		Fourier Map		Distances and Angles	
Quantification			Refinement Control		
Pattern List		Scan List	Peak List		Anchor Scan Data
No.	Pos. [*2Th.]	d-spacing [Å]	FWHM [...]	h	
1	41.6243	2.51940	0.4723		
2	52.8959	2.00984	0.2165		
3	67.8739	1.60341	0.2362		
4	77.5586	1.42819	0.5760		


- d. Dari data XRF diketahui baja yang dianalisis unsur utamanya adalah Fe sehingga pencarian yang dilakukan hanya pada elemen Fe saja. Berikut langkah-langkahnya,



Klik menu Reference Patterns → Retrieve Pattern by → Restrictions..., pencocokan ini akan memakan waktu dan mungkin tidak akan memperoleh hasil.




Klik elemen yang diinginkan, bila di klik lebih dari satu kali warna pilihannya akan berubah sesuai urutan LEGEND



Menunggu.....

ditampilkan entries pada "Pattern List"



Structure Plot | Fourier Map | Distances and Angles

Quantification | Refinement Control

Pattern List | Scan List | Peak List | Anchor Scan Data

Accepted Pattern:

No.	Visible	Ref. Code	Compound N...	Chemical
1	<input checked="" type="checkbox"/>	00-001-1252	Iron	Fe
2	<input checked="" type="checkbox"/>	00-001-1262	Iron	Fe
3	<input checked="" type="checkbox"/>	00-001-1267	Iron	Fe
4	<input checked="" type="checkbox"/>	00-003-1050	Iron	Fe
5	<input checked="" type="checkbox"/>	00-006-0696	ferrite	Fe
6	<input checked="" type="checkbox"/>	00-034-0529	Hexaferrum	Fe
7	<input checked="" type="checkbox"/>	00-050-1275	Iron	Fe
8	<input checked="" type="checkbox"/>	00-052-0513	austenite	Fe
9	<input type="checkbox"/>	01-085-1410	Iron	Fe

Candidates:

No.	Ref. Code	S.	Compound Name	Chemical
-----	-----------	----	---------------	----------

Berikutnya menyisihkan pola difraksi yang tidak sesuai. Dari gambar di atas, yang cocok adalah pola difraksi 00-001-1267.

Peak List

Ref. Code	Color
00-001-1267	Green

Position [ $^{\circ}2\theta$ ]

Structure Plot | Fourier Map | Distances and Angles

Quantification | Refinement Control

Pattern List | Scan List | Peak List | Anchor Scan Data

Accepted Ref. Pattern: 00-006-0696

No.	Visible	Ref. Code	Compound N...	Chemical
1	<input type="checkbox"/>	00-001-1252	Iron	Fe
2	<input type="checkbox"/>	00-001-1262	Iron	Fe
3	<input checked="" type="checkbox"/>	00-001-1267	Iron	Fe
4	<input type="checkbox"/>	00-003-1050	Iron	Fe
5	<input type="checkbox"/>	00-006-0696	ferrite	Fe
6	<input type="checkbox"/>	00-034-0...		
7	<input type="checkbox"/>	00-050-1...		
8	<input type="checkbox"/>	00-052-0...		
9	<input type="checkbox"/>	01-085-1...		

Candidates:

No.	Ref. Code	S.	Compound Name	Chemical
-----	-----------	----	---------------	----------

Remove Pattern

Remove All Patterns

Find...

Toggle visible

Modify RIR...

Simulate Scan from Pattern

Connect Pattern to Structure

Convert Pattern to Phase

Save Pattern List as Subset...

Load Subset to Pattern List...

Autosize Columns

Customize Pattern List...


Copy List

Print List...

Save List as...

Pada Pattern List, hilangkan tanda check list pada kolom Visible untuk pola difraksi yang tidak cocok, kemudian tekan CTRL+A untuk memilih semua pola difraksi, tahan tombol CTRL dan klik kiri pada pola difraksi 00-001-1267 untuk membatalkan pemilihan pada pola difraksi tersebut.

Selanjutnya klik kanan pada pola difraksi yang terpilih, klik pada pilihan Remove Pattern.



e. Memulai Analisis Rietveld

Bila database PDF-4 yang digunakan, pola difraksi pada Pattern List dapat langsung dijadikan acuan untuk analisis Rietveld.

Oleh karena itu diperlukan data kristalografi dari database COD yang sudah berisi data posisi atom. Untuk keperluan analisis diperoleh data CIF dari COD (bisa melihat penggunaan aplikasi MATCH!) dengan nomor kartu 96-900-8537 dan 96-900-6604.

X'Pert HighScore Plus - [bajaST37]

File Edit View Treatment Reference P

Automatic Mode 0 Parameter(s) varied 0 Constraint(s)

Structure Plot Fourier Map Distances and Angles  
 Pattern List Scan List Peak List Anchor Scan Data  
 Quantification Refinement Control

Name	Refine	Value
Global Parameters		

**Selected object: Global Parameters**

**Background**

Method Polynomial

Use Extended Background Terms

Flat Background 0.000000

**Coefficients**

Coefficient 1	0.000000
Coefficient 2	0.000000
Coefficient 3	0.000000
Coefficient 4	0.000000
Coefficient 5	0.000000

**Agreement Indices**

R expected	0.00000
R profile	0.00000
Weighted R profile	0.00000

Displacement Constraints

Flat Plate Geometry

Displacement Constraints  
 Displacement (mm) 0.00000

Coefficient 1.00000

Coefficient [1/cm] 0.00000

Option Correction

File Edit View Treatment Reference P

New Ctrl+N

Open... Ctrl+O

Insert... Ctrl+I

Close

Save Document Ctrl+S

Save As...

Graphic Page Setup...

Print Graphics... Ctrl+P

Print Preview

Graphic Print Setup...

Send To...

Properties...

Exit

**PANalytical X'Pert HighScore Plus**

Insert data into document

Look in: Analysis

Name	Date modified	Type
bajaST37	1/14/2011 7:32 AM	PANalytic
Entry_96-900-6604_Fe	1/14/2011 7:32 AM	Crystallog
Entry_96-900-8537_Fe	6/18/2012 9:45 PM	Crystallog

File name: Entry\_96-900-6604\_Fe

Files of type: All files (\*.\*)

Open as read-only

Open Cancel

X'Pert HighScore Plus - [bajaST37]


File Edit View Treatment Reference P

Automatic Mode 0 Param

Structure Plot Fourier Map Distances and Angles  
 Pattern List Scan List Peak List Anchor Scan Data  
 Quantification Refinement Control



Name	Refine	Value
Global Parameters		
9006603	<input checked="" type="checkbox"/>	

f. Selanjutnya pertama kali "Refine"



Mode refinement:

1. Automatic
2. Semi Automatic
3. Manual Mode



**Selected object: Global Parameters**

**Background**

Method: Polynomial

Use Extended Background Terms:

Flat Background: 263.1420000

**Coefficients**

Coefficient 1	-1.3945600
Coefficient 2	0.0131522
Coefficient 3	0.0000000
Coefficient 4	0.0000000
Coefficient 5	0.0000000

**Agreement Indices**

R expected	6.24960
R profile	6.02153
Weighted R profile	7.68484
D-statistics	0.77957
Weighted D-statistics	1.32052
Goodness of Fit	1.51205

**General Properties**

Sample Geometry: Flat Plate Geometry

Automatic Cell Constraints:

Automatic Anisotropic Displacement Constraints:

Specimen Displacement [mm]: 0.00000

Zero Shift [2Theta]: 0.64097

Polarisation Correction Coefficient: 1.00000

Linear Absorption Coefficient [1/cm]: 0.00000

Use Brindley Microabsorption Correction:

**Refinement**


**U standard**

**V standard**

**W standard**

Nilai Goodness of Fit yang diperoleh 1.51205

Untuk melihat parameter apa saja yang di “Refine” klik kanan pada area “Global Parameter” lalu klik “Show Refined Values/Constraints.....”




No.	Description	Value	Deviation	Last Shift	Minimum	Maximum	Code	Use Min./M...
1	Zero Shift [*2Theta]	0.640974	0.009762	0.000000	-1.000000	1.000000	41.000	<input checked="" type="checkbox"/>
2	Flat background	263.142000	2.509062	0.000040	0.000000	0.000000	11.000	<input type="checkbox"/>
3	Coefficient 1	-1.394560	0.099327	-0.000001	0.000000	0.000000	21.000	<input type="checkbox"/>
4	Coefficient 2	0.013152	0.000891	0.000000	0.000000	0.000000	31.000	<input type="checkbox"/>
5	9006603 Scale Factor	0.017646	0.000419	0.000000	0.000000	2.147484E9	51.000	<input checked="" type="checkbox"/>
6	9006603 Overall B	0.000000	0.180367	-0.300093	0.000000	10.000000	71.000	<input checked="" type="checkbox"/>
7	9006603 Cell a [Å]	2.877747	0.000482	0.000000	0.000000	0.000000	81.000	<input type="checkbox"/>
8	9006603 Profile W	0.140561	0.002163	0.000000	-1.000000	10.000000	61.000	<input checked="" type="checkbox"/>

g. Proses “Refine” lanjutan.

Untuk melakukan proses “Refine” lanjutan pilih mode refinement menjadi semi automatic. Selanjutnya lakukan penambahan checklist pada parameter global ataupun pada parameter phase.

## 2. MEMBUAT LAPORAN

Pada langkah ini akan membuat laporan dari ditampilkan dari XHP.



Berikut lampiran file yang terbentuk

X'Pert HighScore PANalytical	X'Pert HighScore Report 	date: 10/16/13 for internal use only
---------------------------------	--	--

## Fantastic report on Document1



### Distribution:

Name	Function	0.1	0.5	1.0
N.N.-999	Mastering Group Manager (MGM)	✓		
N.N.-998	Product Ontwikkelt Manager (POM)	✓		
N.N.-1	(AS)			✓
N.N.-2	(DS)			✓
N.N.-3	(BS)			✓
T.P.	Supervisor All Coding			✓

Last saved by Jan Setiawan	Draft	file: bajaST37	rev: 1
© 2013 PANalytical B.V. All Rights Reserved		No. of pages: 1/9	

**Table of Contents:**

**TABLE OF CONTENTS:** ..... 2

**ANCHOR SCAN PARAMETERS: (BOOKMARK1)**..... 3

**GRAPHICS: (BOOKMARK2)**..... 4

**PEAK LIST: (BOOKMARK 3)**..... 5

**IDENTIFIED PATTERNS LIST: (BOOKMARK4)** ..... 6


**PLOT OF IDENTIFIED PHASES: (BOOKMARK 5)**..... 6

**DOCUMENT HISTORY: (BOOKMARK 6)**..... 7

**Anchor Scan Parameters:** (Bookmark1)

Sample Identification	ST370
Comment	Material Science
Comment	PC-APD, Diffraction software
Measurement Date / Time	9/5/2009 1:20:00 PM
Raw Data Origin	PHILIPS-binary (scan) (.RD)
Scan Axis	Gonio
Start Position [°2Th.]	20.0000
End Position [°2Th.]	89.9800
Step Size [°2Th.]	0.0200
Scan Step Time [s]	1.0000
Scan Type	Continuous
Offset [°2Th.]	0.0000
Divergence Slit Type	Fixed
Divergence Slit Size [°]	0.2500
Specimen Length [mm]	10.00
Receiving Slit Size [mm]	0.2000
Measurement Temperature [°C]	0.00
Anode Material	Co
K-Alpha1 [Å]	1.78901
K-Alpha2 [Å]	1.79290
K-Beta [Å]	1.62083
K-A2 / K-A1 Ratio	0.50000
Generator Settings	30 mA, 40 kV
Diffractometer Type	PW3710
Diffractometer Number	1
Goniometer Radius [mm]	173.00
Dist. Focus-Diverg. Slit [mm]	91.00
Incident Beam Monochromator	No
Spinning	No

**Graphics:** (Bookmark2)



**Peak List:** (Bookmark 3)


Pos. [°2Th.]	Height [cts]	FWHM [°2Th.]	d-spacing [Å]	Rel. Int. [%]	Tip width [°2Th.]	Matched by
52.7961	2537.10	0.3749	2.01191	100.00	0.4499	00-001-1267
52.9180	2279.90	0.3749	2.01197	89.86	0.4499	
77.5174	301.83	0.3749	1.42883	11.90	0.4499	00-001-1267
77.7153	236.54	0.3749	1.42886	9.32	0.4499	

Report

X'Pert HighScore

**Identified Patterns List:** (Bookmark4)

Visible	Ref. Code	Score	Compound Name	Displacement [°2 $\theta$ .]	Scale Factor	Chemical Formula
*	00-001-1267	29	Iron	0.000	0.969	Fe

**Plot of Identified Phases:** (Bookmark 5)

**Document History:** (Bookmark 6)

## Insert Measurement:

- File name - "bajaST37.RD"
- Modification time - "10/16/2013 10:38:43 AM"
- Modification editor - "Jan Setiawan"

## Search Peaks:

- Minimum significance - "2.00"
- Minimum tip width - "0.01"
- Maximum tip width - "1.00"
- Peak base width - "2.00"
- Method - "Minimum 2nd derivative"
- Modification time - "2/20/2001 11:55:18 AM"
- Modification editor - "PANalytical"

## Insert Structure(s):

- File name - "Entry\_96-900-6604\_Fe.cif"
- Modification time - "10/16/2013 11:09:08 AM"
- Modification editor - "Jan Setiawan"

## Delete Selected Phase:

- Modification time - "10/16/2013 11:09:17 AM"
- Modification editor - "Jan Setiawan"

## Insert Structure(s):

- File name - "Entry\_96-900-6604\_Fe.cif"
- Modification time - "10/16/2013 11:09:25 AM"
- Modification editor - "Jan Setiawan"

## Change 9006603 Use Flag:

- Modification time - "10/16/2013 11:09:54 AM"
- Modification editor - "Jan Setiawan"

## Delete All Phases:

- Modification time - "10/16/2013 11:10:01 AM"
- Modification editor - "Jan Setiawan"

## Insert Structure(s):

- File name - "Entry\_96-900-6604\_Fe.cif"
- Modification time - "10/16/2013 11:16:09 AM"
- Modification editor - "Jan Setiawan"

## Rietveld Refinement:

- Step No. 1
- Title - "Scale factor"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Step No. 2
- Title - "Flat background"
- Min. Shift/ESD - "0.1"

---

- Switch off after usage - "False"
- Step No. 3
- Title - "Zero shift"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Refine spec. displ. - "False"
- Step No. 4
- Title - "Lattice parameters"
- Min. Shift/ESD - "0.8"
- Switch off after usage - "False"
- Minimum weight % - "-1"
- Step No. 5
- Title - "More background"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- No. additional parameters - "2"
- Step No. 6
- Title - "W (Halfwidth)"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Minimum weight % - "1.5"
- Step No. 7
- Title - "Site occupancy factor and B isotropic"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Minimum atomic number - "10"
- Use overall B - "True"
- Refine mixed sites - "False"
- Minimum weight % - "-1"
- Modification time - "10/16/2013 11:20:35 AM"
- Modification editor - "Jan Setiawan"

## Rietveld Refinement:

- Step No. 1
- Title - "Scale factor"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Step No. 2
- Title - "Flat background"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Step No. 3
- Title - "Zero shift"
- Min. Shift/ESD - "0.1"
- Switch off after usage - "False"
- Refine spec. displ. - "False"
- Step No. 4
- Title - "Lattice parameters"
- Min. Shift/ESD - "0.8"
- Switch off after usage - "False"
- Minimum weight % - "-1"
- Step No. 5
- Title - "More background"

---

**Report****X'Pert HighScore**

- Min. Shift/ESD = "0.1"
- Switch off after usage = "False"
- No. additional parameters = "2"
- Step No. 6
- Title = "W (Halfwidth)"
- Min. Shift/ESD = "0.1"
- Switch off after usage = "False"
- Minimum weight % = "1.5"
- Step No. 7
- Title = "Site occupancy factor and B isotropic"
- Min. Shift/ESD = "0.1"
- Switch off after usage = "False"
- Minimum atomic number = "10"
- Use overall B = "True"
- Refine mixed sites = "False"
- Minimum weight % = "-1"
- Modification time = "10/16/2013 11:45:02 AM"
- Modification editor = "Jan Setiawan"

Untuk mengubah penampilan Report sesuai dengan kebutuhan,

